

# Jazzy<sup>®</sup>

## SELECT 6

The Jazzy<sup>®</sup> Select<sup>™</sup> 6 provides maximum front stability, combined with the advantages of in-line motor technology. Active-Trac<sup>®</sup> Suspension and 15 cm (6") front and rear casters offer superior performance indoors and out. An impressive selection of standard convenience features make it the easiest power chair to use.



**Jazzy<sup>®</sup> Select<sup>™</sup> 6** shown in *Jazzy Red* with gray vinyl high-back seat, PG GC controller and foot platform

## SPECIFICATIONS

<b>Weight capacity</b>	136 kg (300 lbs.)
<b>Maximum speed</b>	Up to 6 km/h (4 mph)
<b>Ground clearance</b>	5 cm (2")
<b>Turning radius</b>	52 cm (20.5")
<b>Overall length</b>	93 cm (36.625") without front riggings
<b>Overall width</b>	60 cm (23.5")
<b>Wheels</b>	Drive: 25 cm (10") solid
<b>Front casters</b>	15 cm (6") solid
<b>Rear casters</b>	15 cm (6") solid
<b>Suspension</b>	Active-Trac®
<b>Drive train</b>	Two-motor in-line, mid-wheel drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	40 amp, PG GC controller
<b>Specialty controls</b>	N/A
<b>Battery charger</b>	Off-board, 3 amp
<b>Per-charge range</b>	Up to 24 km (15 miles)
<b>Battery requirements</b>	(2) 12 volt, deep cycle, 35AH
<b>Battery weight</b>	11 kg (24.5 lbs.) each
<b>Base weight</b>	45 kg (100 lbs.)
<b>Standard seat weight</b>	17 kg (37 lbs.) (medium-back)
<b>Max. seat dimensions</b>	51 cm x 51 cm (20" x 20")
<b>Warranty</b>	2-year limited

## COLOUR CHOICES



Jazzy Red



Jazzy Blue

## FEATURES

- 15 cm (6") casters on front and rear for maximum stability
- In-line motor technology for increased efficiency, torque and range
- Articulating rear caster beam
- Easy access charger port and circuit breaker
- Height and depth-adjustable foot platform
- Built with ease of service in mind

## OPTIONS

- NEW depth-adjustable seat provides fast, simple adjustment for enhanced, personalized comfort\*
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

**Pride**<sup>®</sup>  
 Mobility Products  
 Australia Pty. Ltd.



21 Healey Rd.  
 Dandenong, Victoria 3175  
[www.pridemobility.com](http://www.pridemobility.com)

ACN 088 609 661 ABN 66 088 609 661

Tel: 03 9706 4611  
 Fax: 03 9706 4622

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient